


**1. Part Number:** 675 8050P1000

**2. Photograph:**



### 3. Product Information:

Description:

LiquiThief PP - 1000mm Single use,

Use:

disposable liquid sampler

Material of Construction:

Polypropylene (PP) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- Conforms to REACH regulation 1907/2006
- BSE/TSE free

O-ring Seal – virgin Silicone Rubber

- Conforms to FDA CFR 177.2600

Method of Construction:

Ultrasonic welding

Nominal Length

Fully closed:

991mm

Fully extended:

1818mm

Diameter of Outer Tube:

21mm

Diameter of hole in tip:

1.5mm

Width of Handle:

103mm

Nominal Weight of Sampler:

151g

Moulding Environment:

Class 100,000 Cleanroom

Assembly & Primary Packing Environment: Class 10,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Method of Sterilization:

N/A

Quantity per Box:

20

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This product has not been ATEX tested

## 4. Change Control:

This is a single-use, disposable sampler. Repeated use will lead to failure.

## 5. Country of Origin:

## 6. Repeated Use:

## 7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Jul-24