

1. Part Number: 675 8050P500

2. Photograph:



3. Product Information:

Description:	LiquiThief PP - 500mm Single use, disposable liquid sampler
Use:	
Material of Construction:	Polypropylene (PP) - virgin <ul style="list-style-type: none"> • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • Conforms to REACH regulation 1907/2006 • BSE/TSE free <p>O-ring Seal – virgin Silicone Rubber</p> <ul style="list-style-type: none"> • Conforms to FDA CFR 177.2600
Method of Construction:	Ultrasonic welding
Nominal Length	
Fully closed:	523mm
Fully extended:	881mm
Diameter of Outer Tube:	21mm
Diameter of hole in tip:	1.5mm
Width of Handle:	103mm
Nominal Weight of Sampler:	79g
Moulding Environment:	Class 100,000 Cleanroom
Assembly & Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilization:	N/A
Quantity per Box:	20
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested

4. Change Control:

This is a single-use, disposable sampler. Repeated use will lead to failure.

5. Country of Origin:

6. Repeated Use:

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 10-Jul-24