


**1. Part Number:** 675 8051H450

**2. Photograph:**



**3. Product Information:**

Description:

SteriWare Jumbo Pipette 450mm Single

Use:

use, disposable liquid sampler High

Material of Construction:

Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

A blue dye is added to the Tip. The blue dye certified as being FDA acceptable and does not contain any materials of animal origin.

Method of Construction:

Ultrasonic welding

Nominal Length

420mm

Diameter of Outer Tube:

21mm

Diameter of hole in tip:

Handle End (Blue) = 1.5mm  
Product End (Natural) = 3mm

Nominal Weight of Sampler:

50g

Moulding Environment:

Class 100,000 Cleanroom

Assembly & Primary Packing Environment: Class 10,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Method of Sterilisation:

N/A

Number of Samplers per Box:

50

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This product has not been ATEX tested