



1. Part Number: 675 8051H900

2. Photograph:



3. Product Information:

Description: SteriWare Jumbo Pipette 900mm Single

Use: use, disposable liquid sampler High

Material of Construction: Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

A blue dye is added to one of the tips. The blue dye is certified:

- FDA acceptable for food contact
- Free of materials of animal origin

Method of Construction: Ultrasonic welding

Nominal Length 892mm

Diameter of Outer Tube: 21mm

Diameter of hole in tip: Handle End (Blue) = 1.5mm
Product End (Natural) = 3mm

Nominal Weight of Sampler: 103g

Moulding Environment: Class 100,000 Cleanroom

Assembly & Primary Packing Environment: Class 10,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: N/A

Number of Samplers per Box: 50

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

ATEX: This product has not been ATEX tested