

1. Part Number: 675 8060H501

2. Photograph:



3. Product Information:

Description: SteriWare PowderThief - 500mm – 1ml Tip

Use: Single use, disposable powder sampler

Material of Construction: HDPE (High Density Polyethylene) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction: Ultrasonic welding

Nominal Length: 511mm

Diameter of Outer Tube: 21mm

Nominal Tip Size: 1ml

Width of Handle: 103mm

Length 'L': 5mm

Nominal Weight of Sampler: 93g

Moulding Environment: Class 100,000 Cleanroom

Assembly & Primary Packing Environment: Class 10,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)

Method of Sterilisation: N/A

Number of Samplers per Box: 20

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

ATEX: This device has not been ATEX tested