

1. Part Number: 675 8079A600

2. Product Information:

Description: Sterile Bag Stand-Up 600ml

Use: Single use, disposable sampling bags

Material of Construction: Thermoplastic laminate

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU 10/2011
- BSE/TSE free

Closure Type: Attachment with 2 round wires

Bag Dimensions: 152 x 229mm

Thickness of Bag: 89 microns

Nominal Bag Volume: 600ml

Method of Sterilization: Sterility is provided by irradiation. The sterilization dose ensuring a sterility assurance level (SAL) of 10^{-3} has been established according to ISO 11137.

Packaging: 1000 bags per case (4 x 250)

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years

ATEX: This product has not been ATEX tested

Made and packed in Canada.

3. Country of Origin: This is a single-use, disposable product and should not be re-used.

4. Repeated Use:

5. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Jul-24