

1. Part Numbers: 675 8126H30A

2. Photograph:



3. Product Information:

Description	Wide Mouth HDPE Bottle 30ml Amber
Materials of Construction:	Bottle HDPE (conforms to FDA 177.1520 & USP Class VI) Cap PP (conforms to FDA 177.1520 & USP Class VI)
Manufacturing Environment:	Controlled environment within factory
Nominal Volume:	30ml
Brimful Volume:	34ml
Nominal Height:	63mm
Diameter:	37mm
Internal Neck Diameter:	21mm
Weight:	13g
Pack Size:	72 bottles
Method of Sterilisation:	N/A (these bottles are non-sterile)
BSE/TSE Free:	Yes
Supplied with lid screwed on:	Yes
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Nov-24