


1. Part Number: 675 8360T50

2. Photograph:



3. Product Information:

Description:	AccuPod 1.0ml
Use:	Single use, disposable sampling pod
Material of Construction:	High Density Polyethylene (HDPE) – virgin <ul style="list-style-type: none"> • FDA acceptable CFR 177.1520 • EU Regulations 10/2011 compliant • EC regulation 1935/2004 compliant • BSE/TSE free
Method of Construction:	Body = Single piece, injection moulded Lid = Single piece, injection moulded
Nominal Volume:	5.0ml
Length:	61mm
Width:	13mm
Max. Height:	20mm
Weight:	5g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of AccuPods per Box:	50
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested