

1. Part Number: 675 A3041000

2. Photograph:



3. Product Information:

Description	Economy Scoop 1000ml
Nominal Volume:	1000ml
Material of Construction:	304 stainless steel
Method of Construction:	Spot welded (NB. This method of construction will leave crevices)
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	315mm
Body Length:	200mm
Body Diameter:	120mm
Nominal Weight:	360g
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24