

1. Part Number: 675 A304450

2. Photograph:



3. Product Information:

Description	Economy Scoop 450ml
Nominal Volume:	450ml
Material of Construction:	304 stainless steel
Method of Construction:	Spot welded (NB. This method of construction will leave crevices)
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	255mm
Body Length:	145mm
Body Diameter:	75mm
Nominal Weight:	220g
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24