

1. Part Number: 675 A3161250

2. Photo:



3. Product Information:

Description	PharmaScoop® 1250ml
Nominal Volume:	1250ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	305mm
Body Length:	195mm
Body Diameter:	125mm
Nominal Weight:	540g
Total Surface Area:	$1.24 \times 10^5 \text{ mm}^2$ (approx.)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 18-Dec-24