

1. Part Number: 675 A3161700

2. Photo:



3. Product Information:

Description	PharmaScoop® 1700ml
Nominal Volume:	1700ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	355mm
Body Length:	230mm
Body Diameter:	135mm
Nominal Weight:	700g
Total Surface Area:	1.57x10 ⁵ mm ² (approx.)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24