

**1. Part Number:** 675 A330100

**2. Photograph:**



### 3. Product Information:

Description	Industrial Pallet Knife 100mm
Material:	Stainless Steel - grade 420
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	200mm
Blade Length:	100mm
Max. Blade Width:	20mm
Nominal Weight:	60g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

**Note:** Not recommended for GMP use as there is a small crevice where the blade meets the handle.  
For GMP use see the A123 GMP Pallet Knife range.

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24