

1. Part Number: 675 A333100

2. Photograph:



3. Product Information:

Description:

PTFE Scoop 100ml

Material:

PTFE (Polytetrafluoroethylene), virgin

- FDA acceptable
- USP Class VI compliant

Nominal Volume:

100ml

Body Diameter:

65mm

Body Length:

95mm

Overall Length

195mm

Weight:

190g

Recommended Storage Conditions:

Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24