

1. Part Number: 675 A3441500

2. Photograph:



3. Product Information:

Description	Paddle 1500mm
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	1500mm
Diameter of Handle:	25mm
Blade Length:	300mm
Blade Width:	150mm
Nominal Weight:	2.0kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24