

**1. Part Number:** 675 A3462000

**2. Photograph:**



### 3. Product Information:

Description	Scorpion PharmaScoop 2000ml
Nominal Volume:	2000ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Body Length:	300mm
Body Width:	160mm
Maximum Height (inc. handle):	140mm
Weight:	750g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24