

1. Part Number: 675 A656100

2. Photograph:



3. Product Information:

Description	Scraper 316L stainless steel
Material of Construction:	316L stainless steel
Method of Construction:	Fully welded, the weld is ground and polished
Surface Finish:	Better than 0.5 microns Ra
Max. Blade Width:	110mm
Blade Length:	115mm
Overall Length:	235mm
Nominal Weight:	285g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Storage Conditions:	Dry, ambient conditions

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25