

**1. Part Number:** 675 A426100

**2. Photograph:**



### 3. Product Information:

|                                 |                                                                                         |
|---------------------------------|-----------------------------------------------------------------------------------------|
| Description:                    | Mini-Pot 100ml                                                                          |
| Method of Construction:         | Crevice free                                                                            |
| Material of Construction:       | 316L Stainless Steel                                                                    |
| Nominal Volume:                 | 100ml                                                                                   |
| Brimful Volume:                 | 108ml                                                                                   |
| Overall Height (body & lid):    | 73mm                                                                                    |
| Body Height:                    | 68mm                                                                                    |
| Maximum Body Diameter:          | 50mm                                                                                    |
| Lid Diameter:                   | 52mm                                                                                    |
| Nominal Weight:                 | 250g                                                                                    |
| Finish:                         | Bright Polish (Ra < 0.8 microns)                                                        |
| BSE/TSE statement:              | All polishing compounds used in the manufacture of this product are of vegetable origin |
| Recommended Storage Conditions: | Dry and ambient temperature                                                             |

Note: If you require container with a watertight seal please see our containers type A836

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Jan-25