

1. Part Number: 675 A42625

2. Photograph:



3. Product Information:

Description:	Mini-Pot 25ml
Method of Construction:	Crevice free
Material of Construction:	316L Stainless Steel
Nominal Volume:	25ml
Brimful Volume:	29ml
Overall Height (body & lid):	36mm
Body Height:	33mm
Maximum Body Diameter:	38mm
Lid Diameter:	41mm
Nominal Weight:	110g
Finish:	Bright Polish (Ra < 0.8 microns)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: If you require container with a watertight seal please see our containers type A836

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Jan-25