

1. Part Number: 675 A643100

2. Photograph:



3. Product Information:

Description:	Industrial Scoop
Nominal Scoop Volume:	100ml
Material of Construction:	Stainless Steel
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	205mm
Body Length:	117mm
Max. Body Width:	65mm
Max. Body Height:	45mm
Nominal Weight:	100g
Recommended Storage Conditions:	Dry and ambient temperature

Note: Due to the method of manufacture these scoops may not be suitable for hygienic applications. For hygienic applications we recommend our range of PharmaScoops®.

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 28-Aug-24

Note: All dimensions are nominal. Please contact us if any dimension is critical for your application.