


1. Part Number: 675 A81610

2. Photograph:



3. Product Information:

Description

Tri-Ferrule Bottle 10ml

Materials of Construction:

- Body: 316L stainless steel
- Lid: 316L stainless steel
- Gasket: Silicone (FDA 21 CFR 177.2600 compliant)

Nominal Volume:

10ml

Method of Construction:

Crevice free interior, welds ground & polished

Surface Finish:

Better than 0.5 microns Ra

Dimensions:

A = 50.5 mm
 B = 34 mm
 C = 48 mm
 D = 44 mm
 E = 35 mm
 Max. Width over the Clamp = 95 mm

Nominal Weight (body, lid, clamp & gasket): 580 g

BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin.

Recommended Storage Conditions:

Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25