


1. Part Number: 675 8073H500

2. Photograph:



3. Product Information:

Description:

Slot Sampler - 620mm

Use:

Single use, disposable powder sampler

Material of Construction:

High Density Polyethylene (HDPE) - virgin

- Conforms to FDA 21 CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Ultrasonic welding

Overall Length:

620mm

Diameter of Outer Tube:

26mm


Diameter of Inner Blade:

22mm

Nominal Weight of Sampler:

148g

Slot Positions:



Number of Sampling Slots:

3

Sampling Slot Width:

8mm

Assembly & Primary Packing Environment: Class 100,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Method of Sterilisation:

N/A

Number of Samplers per Box:

20