

1. Part Number: 675 A39610

2. Photograph:



3. Product Information:

Description	Long Handled PharmaScoop® 10ml
Nominal Volume:	10ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	345mm
Handle Length:	300mm
Body Length:	45mm
Body Diameter:	25mm
Nominal Weight:	64g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 3-Jan-25