

**1. Part Number:** 675 A3961000

**2. Photograph:**



### 3. Product Information:

Description	Long Handled PharmaScoop® 1000ml
Nominal Volume:	1000ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	835mm
Handle Length:	600mm
Body Length:	235mm
Body Diameter:	110mm
Nominal Weight:	990g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Jan-25