

1. Part Numbers: 675 A1251000

2. Photograph:



3. Product Information:

Description	PTFE Bottle 1000ml with screw on Cap
Material of Construction:	PTFE
Nominal Volume:	1000ml
Overall Diameter:	100mm
Internal Diameter of Neck:	57mm
Overall Height:	185mm
Nominal Weight:	530g (body & lid)
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Dec-24