

**1. Part Numbers:** 675 A12525

**2. Photograph:**



### 3. Product Information:

Description	PTFE Bottle 25ml with screw on Cap
Material of Construction:	PTFE
Nominal Volume:	25ml
Overall Diameter:	33mm
Internal Diameter of Neck:	18mm
Overall Height:	61mm
Nominal Weight:	42g (body & lid)
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:

A handwritten signature in blue ink, appearing to read 'Kuldeep Flora', written over a white background.

Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Dec-24