

1. Part Numbers: 675 A125250

2. Photograph:



3. Product Information:

Description	PTFE Bottle 250ml with screw on Cap
Material of Construction:	PTFE
Nominal Volume:	250ml
Overall Diameter:	67mm
Internal Diameter of Neck:	25mm
Overall Height:	110mm
Nominal Weight:	195g (body & lid)
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Dec-24