

1. Part Number: 675 A26650

2. Photograph:



3. Product Information:

Description

Tight Head Drum 50L (Closed Top Drums)

Materials of Construction:

- Body: 316L stainless steel
- Lids: 316L stainless steel
- Gaskets: EPDM (FDA 177.2600 compliant)

Nominal Volume:

50 Litres

Brimful Volume:

52 Litres

Method of Construction:

Welds ground and polished

Surface Finish:

Better than 0.5 microns Ra

External Diameter

362mm

Height:

568mm

Lid Sizes:

2" and 3/4"

Drum Wall Thickness:

1.2mm

Graduation Marks:

None

Handles:

None

Weight:

5.6 Kg

BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin.

Recommended Storage Conditions:

Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24