


1. Part Number: 675 A9265

2. Photograph:



3. Product Information:

Description

Toggle Drum 5L

Materials of Construction:

- Body: 316L stainless steel
- Lid: 316L stainless steel
- Toggles: 304 stainless steel
- Handles: 304 stainless steel
- Gasket: Silicone (FDA 177.2600 compliant & BSE/TSE free)

Nominal Volume:

5 Litres

Brimful Volume:

6 Litres

Method of Construction:

Welds ground and polished

Surface Finish:

Better than 0.5 microns Ra

Dimensions

A = 196 mm
 B = 216 mm
 C = 200 mm

Graduation Marks:

None

Handles:

2 side handles welded onto patch plates

Number of Toggles:

4

Weight:

2.7 Kg (body + gasket + lid)

BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin.

Recommended Storage Conditions:

Dry and ambient temperature

Notes:

These containers are primarily designed for storing powders. They are not suitable for the storage or transportation of beer.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 20-Sep-24