

**1. Part Number:** 675 5035T1001

**2. Photo:**



**3. Product Information:**

Description

Liquid Master PTFE – 1000mm

Materials of Construction:

PTFE & FEP (for product contact parts)

- FDA 21 CFR 177.1550
- USP Class VI
- EU regulation 10/2011
- EC regulation 1935/2004

Outer Diameter of Outer Tube:

25mm

Internal Diameter of Outer Tube:

19mm

Max. Diameter of Tip:

32mm

Diameter of hole in end of tip:

10mm

Overall Length:

1000mm

BSE/TSE Free:

Yes

Nominal Weight:

585g

Recommended Storage Conditions:

Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Dec-24