


1. Part Number: 675 5035T601

2. Photo:



3. Product Information:

Description

Liquid Master PTFE – 600mm

Materials of Construction:

PTFE & FEP (for product contact parts)

- FDA 21 CFR 177.1550
- USP Class VI
- EU regulation 10/2011
- EC regulation 1935/2004

Outer Diameter of Outer Tube:

25mm

Internal Diameter of Outer Tube:

19mm

Max. Diameter of Tip:

32mm

Diameter of hole in end of tip:

10mm

Overall Length:

600mm

BSE/TSE Free:

Yes

Nominal Weight:

370g

Recommended Storage Conditions:

Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Dec-24