

1. Part Number: 675 5061A250

2. Photo:



Closed Position

Note: Telescopic Rod is sold separately.



Open Position

Note: Telescopic Rod is sold separately.

3. Product Information:

Description	Snapper 250ml – 316Lss
Materials:	Half Spheres: 316L stainless steel Frame: 316L stainless steel Nuts & Bolts: 316L stainless steel Springs: 304 stainless steel
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	250ml
Max. Width when fully open:	275mm
Length:	195mm
Nominal Weight:	630g
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: Suitable for use with telescopic rods (Type 5061T)

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 5-Dec-24