

**1. Part Number:** 675 A2041000

**2. Photograph:**



### 3. Product Information:

Description	Stainless Jug 1000ml
Material:	304 Stainless Steel
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	1000ml
Overall Height:	157mm
Diameter (at top):	97mm
Max Width (inc. handle & spout):	145mm
Nominal Weight:	480g
Graduation Marks?	Yes – every 100ml
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24

**1. Part Number:** 675 A2042000

**2. Photograph:**




### 3. Product Information:

Description	Stainless Jug 2000ml
Material:	304 Stainless Steel
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	2000ml
Overall Height:	172mm
Diameter (at top):	129mm
Max Width (inc. handle & spout):	190mm
Nominal Weight:	700g
Graduation Marks?	Yes – every 100ml
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24

**1. Part Number:** 675 A2043000

**2. Photograph:**



**3. Product Information:**

Description	Stainless Jug 3000ml
Material:	304 Stainless Steel
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	235mm
Diameter (at top):	120mm
Diameter (at bottom):	148mm
Nominal Weight:	1000g
Graduation Marks?	No
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24

**1. Part Number:** 675 A2045000

**2. Photograph:**



### 3. Product Information:

Description	Stainless Jug 5000ml
Material:	304 Stainless Steel
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	5000ml
Overall Height:	265mm
Diameter (at top):	140mm
Diameter (at bottom):	172mm
Nominal Weight:	1580g
Graduation Marks?	No
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24