

1. Part Number: 675 A2042000

2. Photograph:




3. Product Information:

Description	Stainless Jug 2000ml
Material:	304 Stainless Steel
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	2000ml
Overall Height:	172mm
Diameter (at top):	129mm
Max Width (inc. handle & spout):	190mm
Nominal Weight:	700g
Graduation Marks?	Yes – every 100ml
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24