

1. Part Number: 675 A2043000

2. Photograph:



3. Product Information:

Description	Stainless Jug 3000ml
Material:	304 Stainless Steel
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	235mm
Diameter (at top):	120mm
Diameter (at bottom):	148mm
Nominal Weight:	1000g
Graduation Marks?	No
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24

1. Part Number: 675 A2045000

2. Photograph:



3. Product Information:

Description	Stainless Jug 5000ml
Material:	304 Stainless Steel
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	5000ml
Overall Height:	265mm
Diameter (at top):	140mm
Diameter (at bottom):	172mm
Nominal Weight:	1580g
Graduation Marks?	No
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24