

1. Part Number: 675 A204500

2. Photograph:



3. Product Information:

Description	Stainless Jug 500ml
Material:	304 Stainless Steel
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	500ml
Overall Height:	132mm
Diameter (at top):	75mm
Max Width (inc. handle & spout):	115mm
Nominal Weight:	260g
Graduation Marks?	Yes – every 100ml
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24