


**1. Part Number:** 675 A68615L

**2. Photograph:**



### 3. Product Information:

Description	Lid for 15L Bucket with Spout
Material of Construction:	316L stainless steel
Method of Construction:	Handle spot welded to top of lid Crevice free interior
Surface Finish:	Better than 0.5 microns Ra
Max. width over spout:	380mm
Diameter:	330mm
Nominal Weight:	1kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature
Note:	Bucket with Spout 15L (A686-15) is sold separately to the lid.

### 4. Approval:

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 8-Jan-25