

1. Part Numbers: 675 A68620

2. Photograph:



3. Product Information:

Description	Bucket with Spout 20L
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free interior Lugs for handle are spot welded to exterior of bucket
Surface Finish:	Better than 1.0 microns Ra
Nominal Volume:	20 L
Brimful Volume:	20.0 L
Max. width over spout:	375mm
Internal Diameter at Top:	330mm
Diameter at Base:	200mm
Bucket Height:	387mm
Graduations:	No
Skirt Band:	Yes
Nominal Weight:	4.0kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25