


1. Part Number: 675 A68620L

2. Photograph:



3. Product Information:

Description	Lid for 20L Bucket with Spout
Material of Construction:	316L stainless steel
Method of Construction:	Handle spot welded to top of lid Crevice free interior
Surface Finish:	Better than 0.5 microns Ra
Max. width over spout:	400mm
Diameter:	350mm
Nominal Weight:	1.1kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature
Note:	Bucket with Spout 20L (A686-20) is sold separately to the lid.

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 8-Jan-25