

1. Part Number: 675 A406100

2. Photograph:



3. Product Information:

Description	Straight Sided Jug 100ml
Material of Construction:	316L stainless steel
Nominal Volume:	100ml
Brimful Volume:	115ml
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug.
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	70mm
Top Diameter	53mm
Top Diameter – including spout:	61mm
Max. Width – from spout to handle	86mm
Nominal Weight:	145g
Graduation Marks?	None
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Aug-24