

1. Part Number: 675 A4061000

2. Photograph:



3. Product Information:

Description	Straight Sided Jug 1000ml
Material of Construction:	316L stainless steel
Nominal Volume:	1000ml
Brimful Volume:	1200ml
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug.
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	161mm
Top Diameter	105mm
Top Diameter – including spout:	113mm
Max. Width – from spout to handle	152mm
Nominal Weight:	563g
Graduation Marks?	Yes (every 100ml starting from 300ml)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Apr-25

1. Part Number: 675 A4062000

2. Photograph:



3. Product Information:

Description	Straight Sided Jug 2000ml
Material of Construction:	316L stainless steel
Nominal Volume:	2000ml
Brimful Volume:	2600ml
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug.
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	171mm
Top Diameter	150mm
Top Diameter – including spout:	161mm
Max. Width – from spout to handle	194mm
Nominal Weight:	844g
Graduation Marks?	Yes (every 500ml starting from 1000ml)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Aug-24

1. Part Number: 675 A4063000

2. Photograph:



3. Product Information:

Description	Straight Sided Jug 3000ml
Material of Construction:	316L stainless steel
Nominal Volume:	3000ml
Brimful Volume:	3350ml
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug.
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	223mm
Top Diameter	147mm
Top Diameter – including spout:	161mm
Max. Width – from spout to handle	194mm
Nominal Weight:	1080g
Graduation Marks?	Yes (every 500ml starting from 1000ml)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Aug-24

1. Part Number: 675 A4065000

2. Photograph:



3. Product Information:

Description	Straight Sided Jug 5000ml
Material of Construction:	316L stainless steel
Nominal Volume:	5000ml
Brimful Volume:	5470ml
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug.
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	245mm
Top Diameter	178mm
Top Diameter – including spout:	192mm
Max. Width – from spout to handle	225mm
Nominal Weight:	1760g
Graduation Marks?	Yes (every 500ml starting from 1000ml)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Aug-24