

**1. Part Number:** 675 A4063000

**2. Photograph:**



### 3. Product Information:

|                                   |                                                                                         |
|-----------------------------------|-----------------------------------------------------------------------------------------|
| Description                       | Straight Sided Jug 3000ml                                                               |
| Material of Construction:         | 316L stainless steel                                                                    |
| Nominal Volume:                   | 3000ml                                                                                  |
| Brimful Volume:                   | 3350ml                                                                                  |
| Method of Construction:           | Crevice free interior. Handle fully welded onto body of jug.                            |
| Surface Finish:                   | Better than 0.5 microns Ra                                                              |
| Overall Height:                   | 223mm                                                                                   |
| Top Diameter                      | 147mm                                                                                   |
| Top Diameter – including spout:   | 161mm                                                                                   |
| Max. Width – from spout to handle | 194mm                                                                                   |
| Nominal Weight:                   | 1080g                                                                                   |
| Graduation Marks?                 | Yes (every 500ml starting from 1000ml)                                                  |
| BSE/TSE statement:                | All polishing compounds used in the manufacture of this product are of vegetable origin |
| Recommended Storage Conditions:   | Dry and ambient temperature                                                             |

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Aug-24