

**1. Part Number:** 675 A406500

**2. Photograph:**



### 3. Product Information:

Description	Straight Sided Jug 500ml
Material of Construction:	316L stainless steel
Nominal Volume:	500ml
Brimful Volume:	600ml
Method of Construction:	Crevice free interior. Handle fully welded onto body of jug.
Surface Finish:	Better than 0.5 microns Ra
Overall Height:	122mm
Top Diameter	87mm
Top Diameter – including spout:	100mm
Max. Width – from spout to handle	132mm
Nominal Weight:	367g
Graduation Marks?	Yes (every 100ml starting from 200ml)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Aug-24