


1. Part Number: 675 A56630

2. Photograph:



A = 119mm
B = 98mm
C = 602mm
D = 401mm
E = 303mm

3. Product Information:

Description	TC Container 30L
Material of Construction:	Body: 316L stainless steel Side Handles: 304 stainless steel Lid: 316L stainless steel Clamp: 316L stainless steel (may be 304 stainless steel) Gasket: Silicone (conforms to FDA CFR 177.2600)
Method of Construction:	Fully welded, the weld is ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	30L
Handles:	Two side handles welded on to patch plates
Graduation Marks:	None
Pressure Rating:	None
Nominal Weight (body, lid, clamp & gasket):	6.8 Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this container are of vegetable origin
Storage Conditions:	Dry, ambient conditions
Notes:	The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3 These containers are primarily designed for storing powders. They are not suitable for the storage or transportation of beer.

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 25-Sep-24