

1. Part Number: 675 A3661

2. Photograph:



3. Product Information:

Description	Square PharmaScoop
Nominal Volume:	1500ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	315 mm
Body Length:	190 mm
Body Width:	160 mm
Body Height:	48 mm
Nominal Weight:	490 g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 20-Dec-24