

1. Part Number: 675 A5042

2. Photograph:



3. Product Information:

Description	Stainless Spoon 10ml 210mm
Material of Construction:	304 stainless steel
Method of Construction:	Fully welded, the weld is ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	10 ml
Overall Length:	210 mm
Bowl Length:	65 mm
Max. Bowl Width:	32 mm
Weight:	68 g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Storage Conditions:	Dry, ambient conditions

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Jan-25