

1. Part Number: 675 A5043

2. Photograph:



3. Product Information:

Description	Stainless Spoon 2ml 300mm
Material of Construction:	304 stainless steel
Method of Construction:	Fully welded, the weld is ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	2 ml
Overall Length:	300 mm
Bowl Length:	36 mm
Max. Bowl Width:	19 mm
Weight:	45 g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Storage Conditions:	Dry, ambient conditions

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Jan-25