


1. Part Number: 675 8115N150

2. Photograph:



3. Product Information:

Description:	SteriWare V-Spatula Type 2 - 150mm
Use:	Single use, powder spatula
Material of Construction:	Polystyrene - virgin • Conforms to FDA CFR 177.1640 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Blade Shape:	V
Blade Style:	Closed Tip
Overall Length	248mm
Blade Length:	147mm
Blade Width:	10mm
Nominal Weight:	11g
Nominal Volume (flat-filled):	2.3ml
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	N/A
Number of Samplers per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested